



TEOXANE

FOR LIPS
THAT
LOVE TO
MOVE

RHA Kiss Collection
BY TEOXANE



INTRODUCING

RHA® KISS VOLUME* A UNIQUE BALANCE OF **STRETCH, STRENGTH, AND STRUCTURAL INTEGRITY**¹


to support natural movement while increasing volume²



DEGREE OF MODIFICATION²
3.6%
EFFECT DURATION³
12 to 18 months
VOLUME³
2 X 1ml



**RHA KISS VOLUME (PRESERVED
NETWORK TECHNOLOGY: PNT) AND
COMPETITORS RHEOLOGICAL PROPERTIES¹**



Adapts faster and more extensively to movement²

- ● ● ● STRETCH¹
- ● ● ● STRENGTH¹
- ● ● ● STRUCTURAL INTEGRITY¹
(resistance to fragmentation)


Featuring an Extra Thin Wall (ETW) 30G needle for enhanced precision and a gentler injection experience.⁴

*RHA® Kiss Volume Lido has the same formulation and rheology as Teosyal RHA3, but displays an ETW 30G needle in the pack
**OBT™: Optimal Balance Technology

Beautify, rejuvenate
and restore lips

PROVEN EFFICIENCY FEWER TOUCH-UPS, LESS PRODUCT, EFFECTIVE IN FULLER LIPS⁵

01 / Fewer touch-ups⁵




➤ Fewer touch-up sessions for RHA Kiss Volume
Assessed on the TLFS (Teoxane Lip Fullness Scale), from grades 1 to 3

02 / Less product⁵

➤ Clinical studies show subjects treated with **RHA Kiss Volume** required less product to achieve optimal cosmetic results compared to the active comparator

03 / Effective in fuller lips⁵




➤ **78.1% vs 65.9%** of patients achieved a **≥1-grade TLFS improvement** with **RHA Kiss Volume** compared to the comparator.

➤ This is particularly reflected in those with fuller lips, with 62.3% achieving at least **≥1-grade TLFS improvement** versus the comparator brand at just 29.4%

Longer lasting results, clinically confirmed⁵

PERCENTAGE OF PATIENTS SHOWING **SIGNIFICANTLY IMPROVED TO MUCH IMPROVED LIP APPEARANCE** ON THE GAIS** AFTER 1 YEAR ARE AS FOLLOWS:



*RHA® Kiss Volume Lido

** Global Aesthetic Improvement Scale

Natural look
and feel,
maintained
over time⁵

↘ **83%** reported patient
satisfaction with
RHA Kiss Volume
lip treatment vs. 70% with active
comparator technology⁺

↘ **93%** of patients
stated that their
smile still felt
natural⁺⁺

↘ **99%** of patients
still enjoy a
natural feel over
time⁺⁺⁺

+ Patients' satisfaction was measured over
time until 12 months, for the Teoxane Lip
Fullness Scale (TLFS) grades 1 to 3

++ Naturalness of treatment. Natural look of
smile over time, as assessed by subjects

+++ Patients rating a grade ≥ 7/10 on the
natural scale, from 3 to 12 months

RHA Kiss Volume* benefits vs competitors^{4,6,7,8}

	RHA KISS VOLUME	RHA KISS VOLUME CHARACTERISTICS**	COMPETITOR 1	COMPETITOR 2	COMPETITOR 3
Technology	Preserved Network Technology ^{9,10}	Heatless crosslinking in gentle conditions ⁹	OBT ^{TM11}	Vycross ^{TM12}	Hylacross ^{TM12}
Indications	Lips volume, philtrum, marionette lines, oral commissures ³	The whole perioral region ³	Lips ⁶	Lips ⁷	Lips ⁸
Degree of Modification (MoD)	~3,6% ²	Less crosslinked, closer to natural HA and its properties ²	6,6% ¹	5,6% ¹	6,9% ¹
Duration	12 to 18 months ³	Lasts for a longer duration ³	Up to 12 months ¹³	Up to 12 months ¹⁴	Up to 12 months ¹⁵
Stretch ¹	●●○○	Malleable gel adapted to facial dynamism, providing more natural-looking results ²	●○○○	●○○○	●●○○
Strength ¹	●●●○	Higher volumising capacity in dynamic conditions ²	●●○○	●○○○	●●○○
Soluble HA Molecular Weight ¹⁶	Higher than competitors 1, 2 & 3**	Better preserves long HA chains ^{**2} and with +4MDa Free HA addition (the closest to natural HA found in young skin) ⁹	Lower***	Lower***	Lower***
Ergonomy at injection	New rod & plunger, shorter syringe with low extrusion force ¹⁰	Shortest syringes with highest precision ¹⁷	Extra long syringe	Extra long syringe	Extra long syringe








*RHA® Kiss Volume Lido

Data on file, RDRE 2207, 2134, 2016. Mechanical properties of
products. Internal data, R&D measurements.

** Compared to listed competitors

*** Compared to RHA Kiss Volume

TAILORED
SOLUTIONS
FOR LIPS
THAT LOVE
TO MOVE

				
	 RHA 1	 RHA Kiss	 RHA Kiss Volume	[3D] Lip
Indication	Hydration and correction of perioral lines	Subtle lip correction and contouring	Advanced dynamic volumisation, for fuller natural-looking lips	Dermocosmetics Nourishing and plumping lip care ²⁰
Hyaluronic Acid concentration	15mg/ml	23mg/ml	23mg/ml	Key active ingredients: hyaluronic acid and collagen microsphere, shea butter and lipopeptide ²¹
Degree of modification / BDDE crosslinking	1.9%	3.1%	3.6%	N/A
Injection depth	Dermis, Subcutaneous	Dermis, Subcutaneous, Subdermis or Submucosa of the red lip	Dermis, Subcutaneous, Subdermis or Submucosa of the red lip	N/A Dermocosmetic product are not designed to be injected
Needle and cannula	Needle 30G or cannula TSK 25G/50mm	Needle 30G or cannula TSK 25G/50mm	Needle Extra Thin Wall 30G 1/2 or cannula TSK 25G/50mm	N/A Dermocosmetic product are not designed to be injected
Packaging	2 x 1ml	2 x 0.7ml	2 x 1ml	1 x 10ml
Effect duration	12 months minimum	12 to 18 months	12 to 18 months	N/A

*The RHA Kiss Collection is made of Teosyal RHA®1, Teosyal RHA® Kiss, RHA® Kiss Volume Lido and [3D] Lip

Every lip is unique. The **RHA Kiss Collection*** offers targeted products to help patients reach different aesthetic goals with precision and confidence.^{3,18,19,20}

Reach your patients' aesthetic goals with our RHA Kiss Collection*

- A AGING ANATOMY ASSESSMENT
- T TECHNIQUES
- P PRODUCTS

RHA 1



- A Hydration and revitalisation of the cutaneous lip and lip contour
- T 30G^{1/2} Needle, bolus & 50mm 25G cannula, retrograde linear threading with **RHA 1** (dermis)
- P **RHA 1** for lip hydration: **1,2mL**

RHA KISS



- A Improving the contour and correcting symmetry in lips
- T 30G^{1/2} Needle, vertical tenting & bolus with **RHA KISS** (labial subdermis & dermis)
- P **RHA KISS** for lip fullness and contour: **0,60mL**

RHA KISS VOLUME



- A Improving the contour and projection in lips
- T 30G^{1/2} ETW Needle, retrograde linear threading with **RHA KISS VOLUME** (labial subdermis)
- P **RHA KISS VOLUME** for lip fullness: **0,90mL**

RHA KISS VOLUME



- A Improving the contour and projection in fuller lips
- T 30G^{1/2} ETW Needle, retrograde linear threading with **RHA KISS VOLUME** (labial subdermis)
- P **RHA KISS VOLUME** for lip fullness and contour: **0,50mL**

*The RHA Kiss Collection is made of Teosyal RHA®1, Teosyal RHA® Kiss, RHA® Kiss Volume Lido and [3D] Lip



TEOXANE

1. RDRE 2207, 2134, 2016. Mechanical properties of products. Internal data, R&D measurements. **2.** Faivre J. et al. Dermatologic Surgery (2021);47(5): e159-e167. **3.** RHA Kiss Volume - instructions for use. **4.** Pavicic T. et al. Influence of needle size and injection angle on the distribution pattern of facial soft tissue fillers. Journal of cosmetic dermatology 18.5 (2019): 1230-1236. **5.** Weinkle S. et al. Aesthetic surgery journal, sjaf135. 14 Jul. 2025. doi:10.1093/asj/sjaf135. **6.** IFU Restylane Kysse. **7.** IFU Juvederm Volbella. **8.** IFU Juvederm Ultra 3. **9.** Data on file. Manufacturing process of RHA gels, TSA-INS-0146_012. **10.** Vantou et al. 2024 Unravelling the key characteristics of the Preserved Network Technology (PNT): Manufacturing the dynamic RHA fillers Poster. **11.** Restylane brochure - Two complementary technologies. **12.** Curcio, Natalie. HA gel technology for currently available HA fillers. American Academy of Dermatology, p.1. 2020. **13.** Hilton S et al. Dermatol Surg. 2018;44:261-269. **14.** Juvederm Volbella XC - Instruction for Use. **15.** Pinsky, Mark A et al. Juvéderm injectable gel: a multicenter, double-blind, randomized study of safety and effectiveness. Aesthetic surgery journal vol. 28,1 (2008): 17-23. **16.** Data on file. RDRE 2022. Analysis of the soluble HA fraction extracted from gels by HPLC-SEC-MALS-RI 2020 and Internal data - R&D measurement. **17.** Data on file. Teoxane Plunger & Finger Grip Usability Evaluation—Report. March 2022. **18.** IFU RHA®1 = Teosyal RHA®1. **19.** IFU RHA®Kiss = Teosyal RHA®Kiss. **20.** in-vivo test carried out by an independent laboratory on 22 volunteers for 28 days. **21.** In-vivo test on ingredients.

Medical devices for professional use only. Product indications and availability vary from country to country. Please consult product instructions locally approved for more information. Medical devices are regulated health products bearing the CE mark, CE2797, under the Medical Device Regulation (EU) 2017/745.

Contact for materiovigilance notification: medical@teoxane.com. GL-Filler+Dermo-2500021

**Teoxane cosmetics comply with the requirements of the European regulation 1223/2009. Cosmetic products are not designed to be injected.*